

Technical Specification



Tub and Shower Chair

DN7140

Stock Numbers:

- **DN7140 – Tub and Shower Chair**
Adjustable seat height from 17" to 21"

Dimensions:





Materials:

Seat and back rest are made from high density polyethylene.

Feet are made from thermal plastic elastomer.

WARNING:

TO PREVENT SERIOUS OR FATAL INJURY: Ensure product is assembled, adjusted and used properly. FOR USE ONLY BY PERSONS WEIGHING LESS THAN 250 LBS. Failure to follow correct installation procedures can cause the legs to become loose. Check after installation to make sure legs are secure, and check periodically to be sure legs remains secure. Contact Moen immediately at 1-800-882-0116 if any components are missing or damaged. Do not modify or substitute any components. FOR QUESTIONS OR ASSISTANCE WITH ASSEMBLY CONTACT MOEN AT 1-800-882-0116. Consult a healthcare professional for advice and instructions specific to user's needs before use.



Stock Number: _____ Tub and Shower Chair
By Moen Incorporated

MOEN®

25300 AL MOEN DRIVE • NORTH OLMSTED, OH 44070
CONSUMER SUPPORT: (800) 289-6636 • FAX: (800) 848-6636
INTERNATIONAL: (440) 962-2000 • FAX: (440) 962-2726
CUSTOMER SERVICE (USA): 800-882-0116 • FAX: (888) 379-2720
WWW.MOEN.COM
© 2018 Moen Incorporated.